
Apera LabSen® Professional pH Electrodes

Apera Instruments LabSen® Professional pH Electrodes are backed by proprietary sensor technologies and components from Switzerland, designed for a wide variety of applications in scientific research and product quality control with more than 30 models.



Glass Membrane

Glass membrane is the most important part of pH electrodes.

- LabSen® pH electrodes are equipped with 4 types of glass membrane to meet needs in various applications: S membrane, H membrane, HF membrane, and PHY membrane.
- LabSen® glass membrane is highly resistant to general impact (totally differentiates from the traditional fragile glass bulb membrane).
- LabSen® glass membrane with different shapes are shown as helow:



Junction

Junction is the electrolyte interface between reference system and the test samples. LabSen® electrode adopts the following junctions:

- Diaphragm the most frequently used junction type, easy to get blocked by protein-containing or suspension solutions.
- Pore without Diaphragm it is used with solid electrolyte, no clogging, and maintenance-free.
- Moveable sleeve easy to clean, suitable for suspension, emulsion, low ion-concentration solution and non-aqueous solution. The infiltration rate of electrolyte is determined by the tightness of the sleeve during installation.
- PTFE a type of Teflon material with multi pores, hard to be contaminated.

Inner Solution

The inner solution of LabSen® electrode is in a distincive dark blue color. With a special gel treatment, inner solution does not flow and will not cause bubble. The electrode can even work upside down.

Reference System

Besides the routine Ag/AgCl reference electrode, LabSen® Electrodes are more likely to adopt the Long-Life reference system and Sliver Ion Trap reference electrode.



Long-Life Reference System

Long-Life reference system is composed of a glass tube, AgCl, and reference silver wire. The top end of the slim glass tube is stuffed with cotton, which will prevent reaction between AgCl and electrolyte when temperature changes. It improves the stability of reference electorde and service life.

LabSen® 211 Routine pH Electrode =



Application

- Suitable for routine lab tests in regular aqueous solutions
- Shock-proof glass membrane **Features**
 - · Long-Life reference

Material	Lead-Free Glass
Temp. Probe	No
Junction	Movable Sleeve
Reference	Long Life
Electrolyte	3M KCl
Length	130mm
Diameter	12mm
Membrane Shape	Cylindrical

LabSen® 221 Precise pH Electrode =



Application • Suitable for low ion-concentration liquid and viscous

Features • Movable sleeve helps generate quick and stable readings in low ion-concentration samples.

Material	Lead-Free Glass
Temp. Probe	Probe
Junction	Ceramic
Reference	Long Life
Electrolyte	3M KCl
Length	150mm
Diameter	12-6mm
Membrane Shape	Cylindrical

LabSen® 241-6 Semi-micro pH Electrode



Range: 0~14 pH Temp.: 0 to 100 °C (32 to 212°F) Connector: BNC/1m cable

Application • Suitable for small volume liquid (≥0.2ml) and test tube samples.

Features • 6mm measuring tip · Long-Life reference system

Material	Lead-Free Glass
Temp. Probe	No
Junction	Ceramic
Reference	Long Life
Electrolyte	3M KCl
Length	130mm
Diameter	12-3mm
Membrane Shape	Cylindrical

LabSen® 241-3 Micro pH Electrode =



Application • Suitable for micro samples (≥30µL), or testing in centrifuge tubes and NMR tubes.

Features • 3mm measuring tip is Φ 3×70mm.

· Long-Life reference system

Material	Lead-Free Glass
Temp. Probe	No
Junction	Ceramic +Single Pore
Reference	Long Life
Electrolyte	Polymer
Length	100mm
Diameter	12-6mm
Membrane Shape	Spear

LabSen® 251 Glass Spear pH Electrode •



Range: 0 to 14 pH Temp.: 0 to 80 °C (32 to 176°F) Connector: BNC/1m cable

Application • Suitable for soft solid and semi-solid food, e.g. cheese, fruit, vegetables, sushi rice, etc.

Features • Spear tip for direct soil pH measurements

· Solid Polymer electrolyte, maintenance-free

Material	PVC
Temp. Probe	No
Junction	Ceramic +Single Pore
Reference	Long Life
Electrolyte	Polymer
Length	90mm
Diameter	15-6mm
Membrane Shape	Spear

LabSen® 551 Spear pH Electrode



· Suitable for solid and semi-solid medium e.g. soil Application

 Spear tip for direct soil pH measurements **Features**

· Solid Polymer electrolyte, maintenance-free

Material	РОМ
Temp. Probe	No
Junction	PTFE
Reference	Long Life
Electrolyte	3M KCl
Length	110mm
Diameter	12mm
Membrane Shape	Flat

LabSen® 371 Flat pH Electrode =



Application

• Suitable for flat surface measurement such as farics, skin, paper, carpet, etc.

· Flat glass membrane, fast response Features

Material	Lead-free Glass
Temp. Probe	No
Junction	Movable Sleeve
Reference	Silver Ion Trap
Electrolyte	3M KCl
Length	130mm
Diameter	12mm
Membrane Shape	Cylindrical

LabSen® 801 Purified Water pH Electrode =



Application

 Suitable for purified water measurement e.g. RO water, distilled water, deionized water etc.

- Features Movable sleeve helps generate quick and stable measurements in low ion-concentration liquid.
 - Silver ion trap reference system

Material	Lead-Free Glass
Temp. Probe	No
Junction	Ceramic
Reference	Silver Ion Trap
Electrolyte	3M KCl
Length	120mm
Diameter	12mm
Membrane Shape	Half Ball

LabSen® 831 HF pH Electrode =



Range: 0 to 11 pH Temp.: 0 to 100°C (32 to 212°F) Connector: BNC/1m Cable

Application •

Suitable for strong acid solutions and solutions containing Hydrofluoric acid (≥3pH).

- Features Special HF glass membrane, HF corrosion resistance.
 - Silver ion trap reference system, preventing the junction from being contaminated by silver sulfide or protein.

Material	Lead-Free Glass
Temp. Probe	No
Junction	Ceramic *3
Reference	Long Life
Electrolyte	Protelyte
Length	120mm
Diameter	12-6mm
Membrane Shape	Conical

LabSen® 851-1 Viscous pH Electrode



Range: 0 to 14 pH Temp.: 0 to 100°C (32 to 212°F) Connector: BNC/1m Cable

Application • Suitable for cosmetics and viscous samples.`

Features • 3 ceramic diaphragm

- Protelyte reference solution
- Silver ion trap referece system

Electrodes for General Purposes

pH, ORP, and Conductivity


201-C Plastic pH Combination Electrode

- Measuring Range: 0 to 14 pH Junction: Ceramic
- Temp. Range: 0 to 80°C
- Reference: Ag/AgCl
- Dimension: Ø12*160mm
- Connector: BNC

Features: Ideal for both lab and in-field use. Gel KCI Electrolyte, no need to refill; Detachable probe cap, easy to clean; Not suitable for testing in strong base solution (pH>12), erosive solutions, or constant testing in high temperature (>60°C)

201T-F Plastic 3-in-1 pH Combination Electrode

In addition to the features of 201-C, it has a built-in thermistor, which allows simultaneous temperature measuring and auto temperature compensation.



301Pt-C Plastic Combination ORP Electrode

- Junction: Ceramic
- Dimension: ø12*160mm
- Sensor: Ø1*5mm Platinum
- Reference: Ag/AgCl
- Connector: BNC

Features: PC housing, Gel KCl electrolyte, no need to refill. Suitable for use in general water solutions and waste water.



2301-C Plastic Conductivity Electrode

- Measuring Range: 0.5 µS/cm to 200 mS/cm
- Electrode Constant: 1.0±0.2 cm⁻¹
- Connector: BNC • Dimension: ø12*155mm
- Sensor: Brush-Resistant Platinum Black Rods

Features: The brush-resistant Platinum black sensor ensures high accuracy in wide measuring ranges. Suitable for lab and field use.

2301T-F Plastic Cond./Temp. Electrode

In addition to the features of 2301-C, it has a built-in thermistor, which allows simultaneous temperature measuring and auto temperature compensation.



2310-C Plastic Conductivity Electrode

- Measuring Range: 20 to 2000 mS/cm
- Electrode Constant: 10± 1 cm⁻¹
- Dimension: Ø12*155mm Connector: BNC
- Sensor: ø5*5 Platinum Black Ring Sensor

Features: High-Concentration Conductivity electrode. Accuracy without calibration: $\leq \pm 10\%$ of readings; Accuracy after calibration: $\leq \pm 1.5\%$ F.S. Suitable for high concentrated electrolyte, sea water, and high concentrated salt water.

2310T-F Plastic Cond. / Temp. Electrode

In addition to the features of 2310-C, it has a built-in thermistor, which allows simultaneous temperature measuring and auto temperature compensation.



2401-C Glass Conductivity Electrode

- Measuring Range: 0.5 μS/cm to 200 mS/cm
- Electrode Constant: 1.0± 0.2 cm⁻¹
- Dimension: ø12*145mm
 Connector: BNC
- Sensor: ø5*5mm Platinum Black

Features: The cavity structure renders higher accuracy and better repeatability, making it suitable for high-precision lab testing.

2401T-F Glass Cond. / Temp. Electrode

In addition to the features of 2401-C, it has a built-in thermistor, which allows simultaneous temperature measuring and auto temperature compensation.



DJS-0.1-C Pure Water Conductivity Electrode

- Measuring Range: 0 to 200 µS/cm
- Electrode Constant: 0.1± 0.02 cm⁻¹
- Dimension: ø12*155mm Connector: BNC
- Sensor: ø7*18mm Platinum Black

Features: equipped with a removable glass flow cell, making it suitable for measurements in pure and ultra-pure water.

DJS-0.1-F Pure Water Cond. / Temp. Electrode

In addition to the features of DJS-0.1-C, it has a built-in thermistor, which allows simultaneous temperature measuring and auto temperature compensation.

