

## 10 μL Gastight Syringe Model 1701 RN, Small Removable Needle, 26s gauge, 2 in, point style 2

Specification Sheet (Part/REF # 80030)



Split/splitless injection consists of a heated chamber with a glass liner into which the sample is injected through the septum. A microsyringe is used to inject the sample through a rubber septum into a flash vaporizer chamber at the head of the column.

## **Product Specifications**

Volume	10 μL
Syringe Type	Gastight
Syringe Series	1700
Termination	Removable Needle (RN)
Needle	Small RN
Gauge	26s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	115°C / 239°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.018 inches (0.460 mm)



Visit <u>www.hamiltoncompany.com/support</u> to find a representative in your area or call us:

US: +1-800-648-5950 EU: +40-356-635-055



Specifications are subject to change without notice		Spec. Version B
Sales Quantity	1 syringe per pack	
Std Needle Dead Volume	0.64 μL	
Needle Outer Diameter	0.0187 inches (0.47 mm)	
Needle Inner Diameter	0.005 inches (0.13 mm)	
Barrel Outer Diameter	0.260 inches (6.60 mm)	

https://www.hamiltoncompany.com/laboratory-products/syringes/80030



US: +1-800-648-5950 EU: +40-356-635-055