## HAMILT®N

## 10 μL, Model 1701 N GelLoading SYR, Cemented NDL,0.3 mm, 1 in, point style 3

Specification Sheet (Part/REF # 84505)



Gel Loading Syringes are designed with small gauge needles that are perfect for loading the lanes on an electrophoresis gel. The most common application is for multiplex DNA sequencing methodologies.

## **Product Specifications**

| -                     |                               |
|-----------------------|-------------------------------|
| Volume                | 10 µL                         |
| Syringe Type          | Gastight                      |
| Syringe Series        | 1700                          |
| Termination           | Cemented Needle (N)           |
| Needle                | Cemented Needle               |
| Gauge                 | 0.3 mm                        |
| Point Style           | 3                             |
| Needle Length         | 1.00 in (25.4 mm)             |
| Fluid Path            | SST, borosilicate glass, PTFE |
| Autoclavable          | No                            |
| Gas Sterilizable      | Yes                           |
| Pressure              | 1000 psig (68.9 bar)          |
| Maximum Temperature   | 50°C / 122°F                  |
| Minimum Temperature   | 10°C / 50°F                   |
| Barrel Inner Diameter | 0.018 inches (0.460 mm)       |
| Barrel Outer Diameter | 0.260 inches (6.60 mm)        |
|                       |                               |



Visit <u>www.hamiltoncompany.com/support</u> to find a representative in your area or call us: US: +1-800-648-5950 EU: +40-356-635-055



| Needle Inner Diameter                               | 0.0025 inches (0.064 mm) |                 |
|---|--------------------------|-----------------|
| Needle Outer Diameter                               | 0.010 inches (0.3 mm)    |                 |
| Std Needle Dead Volume                              | 0.35 μL                  |                 |
| Sales Quantity                                      | 1 syringe per pack       |                 |
| Specifications are subject to change without notice |                          | Spec. Version B |

Specifications are subject to change without notice

https://www.hamiltoncompany.com/laboratory-products/syringes/84505

HAMILT®N

Visit **www.hamiltoncompany.com/support** to find a representative in your area or call us: US: +1-800-648-5950 EU: +40-356-635-055