

## 50 μL Gastight Syringe Model 1805 N, Cemented Needle, 22s gauge, 2 in, point style 2

Specification Sheet (Part/REF # 84881)



Split/splitless injection consists of a heated chamber with a glass liner into which the sample is injected through the septum. A microsyringe is used to inject the sample through a rubber septum into a flash vaporizer chamber at the head of the column.

## **Product Specifications**

Volume	50 μL
Syringe Type	Gastight
Syringe Series	1800
Termination	Cemented Needle (N)
Needle	Cemented Needle
Gauge	22s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.029 inches (0.729 mm)



US: +1-800-648-5950 EU: +40-356-635-055



Barrel Outer Diameter	0.305 inches (7.75 mm)	
Needle Inner Diameter	0.0059 inches (0.15 mm)	
Needle Outer Diameter	0.0283 inches (0.72 mm)	
Std Needle Dead Volume	1.13 μL	
Sales Quantity	1 syringe per pack	
Specifications are subject to change without notice		Spec. Version B

https://www.hamiltoncompany.com/laboratory-products/syringes/84881



US: +1-800-648-5950 EU: +40-356-635-055