

Cat. No	Description	Articles
---------	-------------	----------

<Continued from... **DAIHAN ThermoStable™ IS-30** Precise Shaking Incubators >

Optional Platform/Spring Rack/Flask Holders

- DH.WIS504630 (1) Universal Platform, w485×h405mm, “SP630”, for “IS-30(R)”
- DH.WIS504530 (2) Universal Spring Rack w498×d408×h135mm, “SR530”, for “IS-30(R)”
- DH.WIS504030 (3) Flask Holder, for 50~100ml Flask, “FH100”, for “SHO”, “SHR” and “IS”
- DH.WIS504035 (4) Flask Holder, for 200~300ml Flask, “FH250”, for “SHO”, “SHR” and “IS”
- DH.WIS504040 (5) Flask Holder, for 500~1000ml Flask, “FH500”, for “SHO”, “SHR” and “IS”

Validation (IQ, OQ) Service

- DH.WIS33010 Validation Service(IQ, OQ), “ISV3010”, for “ThermoStable IS-30”
- DH.WIS33011 “ISV3011”, for “ThermoStable IS-30R”



Precise Shaking Incubator, Front Door Type “ThermoStable IS-30”

WiseRemote® Remote Control Software


Connection of Wisd Temperature Devices with Interface RS232C to a PC Export Recorded Data (Excel Spreadsheet)




Main Screen




RS232C Cable & Software



(1) Universal Platform “SP630” (Included)



(2) Universal Spring Rack “SR530” (Optional)



Flask Holders (Optional)